

Hepatitis E Virus ORF3 Antibody FITC Conjugated

Catalog No: #C00152F

Package Size: #C00152F 100ul

Description

Product Name	Hepatitis E Virus ORF3 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	HEV
Crossing Reactivity	HEV
Immunogen Description	KLH conjugated synthetic peptide aa 85-104 114 derived from HEV ORF3
Conjugates	FITC
Target Name	Hepatitis E Virus ORF3
Other Names	structural protein; immunogenic protein; ORF-3; ORF 3; Protein ORF3; pORF3; phosphoprotein; ORF3 [Hepatitis E virus]; ORF3_HEVHY.
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Hepatitis E Virus ORF3

Note: This product is for in vitro research use only and is not intended for use in humans or animals.