

Cytomegalovirus pp65 Antibody HRP Conjugated

Catalog No: #C00195H

Package Size: #C00195H 100ul

Description

Product Name	Cytomegalovirus pp65 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu HCMVPP65
Crossing Reactivity	HCMV PP65
Immunogen Description	KLH conjugated synthetic peptide aa 240-290 561 derived from Cytomegalovirus pp65
Conjugates	HRP
Target Name	Cytomegalovirus pp65
Other Names	65 kDa phosphoprotein; pp65; 65 kDa matrix phosphoprotein; Phosphoprotein UL83; Tegument protein UL83; UL83
Accession No.	Swiss-Prot#Q6SW59NCBI Gene ID3077579
Cell Localization	Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Counteracts the host antiviral immune response by preventing IRF3 from entering the nucleus once activated and phosphorylated. Participates also in the transactivation of viral major immediate-early genes by recruiting host IFI16 to their promoters.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.