

IGBP1 Antibody HRP Conjugated

Catalog No: #C00212H

Package Size: #C00212H 100ul

Description

Product Name	IGBP1 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 156-180 340 derived from human IGBP1
Conjugates	HRP
Target Name	IGBP1
Other Names	IBP1; ALPHA-4; Immunoglobulin-binding protein 1; B-cell signal transduction molecule alpha 4; Protein alpha-4; CD79a-binding protein 1; Protein phosphatase 2 4 6 regulatory subunit; Renal carcinoma antigen NY-REN-16; IGBP1
Accession No.	Swiss-Prot#P78318NCBI Gene ID3476
Cell Localization	Helical
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Associated to surface IgM-receptor; may be involved in signal transduction. Involved in regulation of the catalytic activity of the phosphatases PP2A, PP4 and PP6 by protecting their partially folded catalytic subunits from degradative polyubiquitination until they associate with regulatory subunits.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.