

factor VIII FVIII Antibody FITC Conjugated

Catalog No: #C00223F

Package Size: #C00223F 100ul

Description

Product Name	factor VIII FVIII Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 1470-1520 2351 derived from human factor VIII
Conjugates	FITC
Target Name	factor VIII FVIII
Other Names	AHF; F8B; F8C; HEMA; FVIII; DXS1253E; Coagulation factor VIII; Antihemophilic factor; Procoagulant component; F8
Accession No.	Swiss-Prot#P00451NCBI Gene ID2157
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Factor VIII, along with calcium and phospholipid, acts as a cofactor for factor IXa when it converts factor X to the activated form, factor Xa.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.