

CRF Antibody HRP Conjugated

Catalog No: #C00244H

Package Size: #C00244H 100ul

Description

Product Name	CRF Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 148-196 196 derived from human CRF
Conjugates	HRP
Target Name	CRF
Other Names	CRF; Corticoliberin; Corticotropin-releasing factor; Corticotropin-releasing hormone; CRH
Accession No.	Swiss-Prot#P06850NCBI Gene ID1392
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

This hormone from hypothalamus regulates the release of corticotropin from pituitary gland.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.