

MMP-3 Antibody HRP Conjugated

Catalog No: #C00258H

Package Size: #C00258H 100ul

Description

Product Name	MMP-3 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 430-477 477 derived from human MMP3
Conjugates	HRP
Target Name	MMP-3
Other Names	SL-1; STMY; STR1; CHDS6; MMP-3; STMY1; Stromelysin-1; Matrix metalloproteinase-3; Transin-1; MMP3
Accession No.	Swiss-Prot#P08254NCBI Gene ID4314
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Can degrade fibronectin, laminin, gelatins of type I, III, IV, and V; collagens III, IV, X, and IX, and cartilage proteoglycans. Activates procollagenase.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.