

TFF3 Antibody PE Conjugated

Catalog No: #C00319P

Package Size: #C00319P 100ul

Description

Product Name	TFF3 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 50-80 80 derived from human TFF3
Conjugates	PE
Target Name	TFF3
Other Names	ITF; P1B; TFI; Trefoil factor 3; Intestinal trefoil factor; hITF; Polypeptide P1.B; hP1.B; TFF3
Accession No.	Swiss-Prot#Q07654NCBI Gene ID7033
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.