

IL-13 Antibody HRP Conjugated

Catalog No: #C00334H

Package Size: #C00334H 100ul

Description

Product Name	IL-13 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-75 131 derived from mouse IL-13
Conjugates	HRP
Target Name	IL-13
Other Names	IL-13; Interleukin-13; T-cell activation protein P600; IL13
Accession No.	Swiss-Prot#P20109NCBI Gene ID16163
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Cytokine. Inhibits inflammatory cytokine production. Synergizes with IL2 in regulating interferon-gamma synthesis. May be critical in regulating inflammatory and immune responses (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.