

NF-M Antibody PE Conjugated

Catalog No: #C00405P

Package Size: #C00405P 100ul

Description

Product Name	NF-M Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 120-170 846 derived from rat NF-M
Conjugates	PE
Target Name	NF-M
Other Names	Nfm; Nef3; Neurofilament medium polypeptide; NF-M; 16 kDa neurofilament protein; Neurofilament 3; Neurofilament triplet M protein; Nefm
Accession No.	Swiss-Prot#P12839NCBI Gene ID24588
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Neurofilaments usually contain three intermediate filament proteins: L, M, and H which are involved in the maintenance of neuronal caliber.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.