

CD56 Antibody FITC Conjugated

Catalog No: #C00451F

Package Size: #C00451F 100ul

Description

Product Name	CD56 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 640-690 858 derived from human CD56
Conjugates	FITC
Target Name	CD56
Other Names	CD56; NCAM; MSK39; Neural cell adhesion molecule 1; N-CAM-1; NCAM-1; NCAM1
Accession No.	Swiss-Prot#P13591NCBI Gene ID4684
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

This protein is a cell adhesion molecule involved in neuron-neuron adhesion, neurite fasciculation, outgrowth of neurites, etc.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.