

RRAGC Antibody FITC Conjugated

Catalog No: #C00588F

Package Size: #C00588F 100ul

Description

Product Name	RRAGC Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human RRAGC
Conjugates	FITC
Target Name	RRAGC
Other Names	RAGC; TIB929; YGR163W; GTPase interacting protein 2; GTPase-interacting protein 2; GTR2; Rag C; Rag C protein; Ras related GTP binding C; Ras-related GTP-binding protein C; RP23 29H22.4; RRAGC; RRAGC_HUMAN.
Accession No.	NCBI Gene ID8290
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1ug/Test ICC=1:50-200 IF=1:50-200

Background

This gene encodes a member of the GTR RAG GTP-binding protein family. The encoded protein is a monomeric guanine nucleotide-binding protein which forms a heterodimer with RRAGA and RRAGB and is primarily localized to the cytoplasm. The encoded protein promotes intracellular localization of the mTOR complex. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Dec 2012].

Note: This product is for in vitro research use only and is not intended for use in humans or animals.