

KDM1 Antibody FITC Conjugated

Catalog No: #C00633F

Package Size: #C00633F 100ul

Description

Product Name	KDM1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human KDM1
Conjugates	FITC
Target Name	KDM1
Other Names	Amine oxidase flavin containing domain 2; Amine oxidase flavin containing domain 2; Amine oxidase flavin containing domain protein 2; AOF 2; AOF2; BHC110; BRAF35 HDAC complex protein BHC110; BRAF35-HDAC complex protein BHC110; FAD binding protein BRAF35 HDAC complex, 110 kDa subunit; Flavin contain
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Gene activation and repression is specifically regulated by histone methylation status at distinct lysine residues. Lysine specific demethylase 1 (KDM1 LSD1) is a long-sought histone demethylase that specifically demethylates mono and di methyl histone H3 at K4 and K9. Thus KDM1 is a specific tag for epigenetic transcriptional activation, thereby acting as a corepressor.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.