

CD5 Antibody FITC Conjugated

Catalog No: #C00654F

Package Size: #C00654F 100ul

Description

Product Name	CD5 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 250-300 495 derived from human CD5
Conjugates	FITC
Target Name	CD5
Other Names	T1; LEU1; T-cell surface glycoprotein CD5; Lymphocyte antigen T1 Leu-1; CD5
Accession No.	Swiss-Prot#P06127NCBI Gene ID921
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

May act as a receptor in regulating T-cell proliferation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.