

CD276 B7H3 Antibody PE Conjugated

Catalog No: #C00884P

Package Size: #C00884P 100ul

Description

Product Name	CD276 B7H3 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC,IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 120-170 534 derived from human CD276 B7H3
Conjugates	PE
Target Name	CD276 B7H3
Other Names	B7H3; B7-H3; B7RP-2; 4lg-B7-H3; CD276 antigen; B7 homolog 3; Costimulatory molecule; CD276; PSEC0249; UNQ309 PRO352
Accession No.	Swiss-Prot#Q5ZPR3NCBI Gene ID80381
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

May participate in the regulation of T-cell-mediated immune response. May play a protective role in tumor cells by inhibiting natural-killer mediated cell lysis as well as a role of marker for detection of neuroblastoma cells. May be involved in the development of acute and chronic transplant rejection and in the regulation of lymphocytic activity at mucosal surfaces. Could also play a key role in providing the placenta and fetus with a suitable immunological environment throughout pregnancy. Both isoform 1 and isoform 2 appear to be redundant in their ability to modulate CD4 T-cell responses. Isoform 2 is shown to enhance the induction of cytotoxic T-cells and selectively stimulates interferon gamma production in the presence of T-cell receptor signaling.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.