

Sidekick 1 Antibody

Catalog No: #C00981

Package Size: #C00981-1 100ul

Description

Product Name	Sidekick 1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	affinity purified by Protein A
Applications	IHC;ICC/IF
Species Reactivity	Human; Mouse; Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human Sidekick 1
Accession No.	Swiss-Prot:Q7Z5N4; Gene ID:221935;
Calculated MW	242 kDa
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC 1:100-1:500; ICC/IF 1:100-1:500;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.