

SNAP25 Antibody HRP Conjugated

Catalog No: #C01114H

Package Size: #C01114H 100ul

Description

Product Name	SNAP25 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB,,IHC-P,IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 168-206 206 derived from human SNAP-25
Conjugates	HRP
Target Name	SNAP25
Other Names	RIC4; SEC9; SNAP; RIC-4; SNAP-25; bA416N4.2; dJ168F16.2; Synaptosomal-associated protein 25; Super protein; SUP; Synaptosomal-associated 25 kDa protein; SNAP25
Accession No.	Swiss-Prot#P60880NCBI Gene ID6616
Cell Localization	Cytoplasm, Cell membrane
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

t-SNARE involved in the molecular regulation of neurotransmitter release. May play an important role in the synaptic function of specific neuronal systems. Associates with proteins involved in vesicle docking and membrane fusion. Regulates plasma membrane recycling through its interaction with CENPF.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.