

Somatostatin Antibody HRP Conjugated

Catalog No: #C01124H

Package Size: #C01124H 100ul

Description

Product Name	Somatostatin Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	,IHC-P,IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 89-116 116 derived from human Somatostatin
Conjugates	HRP
Target Name	Somatostatin
Other Names	SMST; Somatostatin; Growth hormone release-inhibiting factor; SST
Accession No.	Swiss-Prot#P61278NCBI Gene ID6750
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IHC-P=1:50-200 IHC-F=1:50-200

Background

Somatostatin inhibits the release of somatotropin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.