

CRBN Antibody PE Conjugated

Catalog No: #C01432P

Package Size: #C01432P 100ul

Description

Product Name	CRBN Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	FCM
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human CRBN
Conjugates	PE
Target Name	CRBN
Other Names	Cereblon; DKFZp781K0715; MGC27358; MRT2A; OTTHUMP00000209555; piL; Protein cereblon; Protein x 0001; 2610203G15Rik; 2900045O07Rik; AF229032; AW108261; CRBN_HUMAN.
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

FCM: 1:20-100

Background

CRBN is a 442 amino acid protein which is highly concentrated in human brain tissue. CRBN functions are thought to be related to energy metabolism, learning and memory. Localized to the cytoplasm, CRBN acts as a protease in mitochondria and is thought to regulate the assembly of KCNT1, as well as the surface expression of KCNT1 in brain regions known to affect memory and learning, such as the hippocampus. The gene encoding CRBN belongs to a family of ATP-dependent ion proteases that play a role in membrane trafficking and proteolysis. Defects in the CRBN gene are associated with mild mental retardation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.