

TNFRSF18 Antibody FITC Conjugated

Catalog No: #C01448F

Package Size: #C01448F 100ul

Description

Product Name	TNFRSF18 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 185-225 228 derived from mouse TNFRSF18
Conjugates	FITC
Target Name	TNFRSF18
Other Names	AITR; Gitr; Tumor necrosis factor receptor superfamily member 18; Glucocorticoid-induced TNFR-related protein; CD357; Tnfrsf18
Accession No.	Swiss-Prot#O35714NCBI Gene ID21936
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Receptor for TNFSF18. Seems to be involved in interactions between activated T-lymphocytes and endothelial cells and in the regulation of T-cell receptor-mediated cell death. Mediated NF-kappa-B activation via the TRAF2 NIK pathway (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.