

KIRREL3 Antibody FITC Conjugated

Catalog No: #C01547F

Package Size: #C01547F 100ul

Description

Product Name	KIRREL3 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human KIRREL3
Conjugates	FITC
Target Name	KIRREL3
Other Names	Kin of IRRE like 3 Drosophila; Kin of IRRE like 3; Kin of irregular chiasm like protein 3 ; KIRRE; KIRREL 3; KIRREL-3; MRD4; NEPH2; Nephrin like 2; PRO19814; PRO4502; yUNQ5923; KIRR3_HUMAN.
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

NEPH2 is a 778 amino acid single-pass type I membrane protein that belongs to the nephrin-like protein family and immunoglobulin superfamily. Expressed in both fetal and adult brain, as well as podocytes of kidney glomeruli, NEPH2 contains five Ig-like C2-type (immunoglobulin-like) domains and is thought to play a role in the hematopoietic supportive capacity of stroma cells. NEPH2 undergoes alternative splicing to produce two isoforms and contains a C-terminal cytoplasmic domain which it uses to interact with Podocin, a podocyte protein involved in ultrafiltration. Defects in the gene encoding NEPH2 are associated with mental retardation autosomal dominant type 4 (MRD4).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.