

SLITRK6 Antibody FITC Conjugated

Catalog No: #C01630F

Package Size: #C01630F 100ul

Description

Product Name	SLITRK6 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human SLITRK6
Conjugates	FITC
Target Name	SLITRK6
Other Names	4832410J21Rik; MGC119595; MGC119596; MGC119597; OTTHUMP00000066012; SLIK6_HUMAN; SLIT and NTRK like family member 6; SLIT and NTRK like protein 6; SLIT and NTRK-like protein 6; Slit and trk like 6; Slit and trk like gene 6; SLITRK 6; Slitrk6.
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

The leucine-rich (LRR) repeat is a 20-30 amino acid motif that forms a hydrophobic α ?? horseshoe fold, allowing it to accommodate several leucine residues within a tightly packed core. All LRR repeats contain a variable segment and a highly conserved segment, the latter of which accounts for 11 or 12 residues of the entire LRR motif. SLITRK6 (SLIT and NTRK-like family, member 6) is an 841 amino acid single-pass type I membrane protein that contains 16 LRR repeats and belongs to the SLITRK family. Expressed at high levels in select brain tissue, as well as in both adult and fetal lung, SLITRK6 functions to suppress neurite outgrowth, playing a role in the regulation of neuronal function. Multiple isoforms of SLITRK6 exist due to alternative splicing events.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.