

GLI1 Antibody FITC Conjugated

Catalog No: #C01706F

Package Size: #C01706F 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GLI1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 610-660 1106 derived from human GLI1
Conjugates	FITC
Target Name	GLI1
Other Names	GLI; Zinc finger protein GLI1; Glioma-associated oncogene; Oncogene GLI; GLI1
Accession No.	Swiss-Prot#P08151NCBI Gene ID2735
Cell Localization	Cytoplasm, Nucleus
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Acts as a transcriptional activator. May regulate the transcription of specific genes during normal development. May play a role in craniofacial development and digital development, as well as development of the central nervous system and gastrointestinal tract. Mediates SHH signaling and thus cell proliferation and differentiation.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.