

CD28 Antibody PE Conjugated

Catalog No: #C02234P

Package Size: #C02234P 100ul

Description

Product Name	CD28 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	FCM
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 20-45 218 derived from mouse CD28
Conjugates	PE
Target Name	CD28
Other Names	T-cell-specific surface glycoprotein CD28; CD28
Accession No.	Swiss-Prot#P31041NCBI Gene ID12487
Cell Localization	Extracellular, Cell membrane
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

FCM: 1:20-100

Background

Involved in T-cell activation, the induction of cell proliferation and cytokine production and promotion of T-cell survival.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.