

DMC1 Antibody FITC Conjugated

Catalog No: #C02255F

Package Size: #C02255F 100ul

Description

Product Name	DMC1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human DMC1
Conjugates	FITC
Target Name	DMC1
Other Names	disrupted meiotic cDNA 1 homolog; dJ199H16.1; DMC 1; DMC1 dosage suppressor of mck1 homolog; DMC1 dosage suppressor of mck1 homolog meiosis specic homologous recombination yeast; DMC1 homologue; DMC1H; HsLim15; LIM15; Meiotic recombination protein DMC1 LIM15 homolog; MGC150472; MGC150473; DMC1_HUMAN
Accession No.	NCBI Gene ID11144
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

DNA repair proteins are necessary for the maintenance of chromosome integrity and are involved in the elimination of premutagenic lesions from DNA. The DNA repair proteins Rad51 and Rad52 are key components of the double-strand-break repair (DSBR) pathway. Rad51 is essential for mitotic and meiotic recombination, and its mutation in yeast and mammalian cells results in chromosome loss. Overexpression of Rad52 confers resistance to ionizing radiation and induces homologous intrachromosomal recombination. Rad52 is thought to be involved in an early stage of Rad51-mediated recombination. Additional proteins involved in the pathway include Nibrin and Dmc1. Nibrin, which complexes with Mre11 and Rad50, is absent in Nijmegen breakage syndrome (NBS) patients. Dmc1 is specifically involved in meiotic recombination. An alternative spliced form of Dmc1, designated Dmc1-D, is deleted for a region between the two motifs involved in nucleotide binding. The alternatively spliced Dmc1-D transcript is detected in both male and female germ cells, indicating that the encoded protein may have a role in mammalian genetic recombination in meiosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.