

MCL1 Antibody FITC Conjugated

Catalog No: #C02650F

Package Size: #C02650F 100ul

Description

Product Name	MCL1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 210-248 350 derived from human MCL1
Conjugates	FITC
Target Name	MCL1
Other Names	TM; EAT; MCL1L; MCL1S; Mcl-1; BCL2L3; MCL1-ES; bcl2-L-3; mcl1 EAT; Induced myeloid leukemia cell differentiation protein Mcl-1; Bcl-2-like protein 3; Bcl-2-related protein EAT mcl1; MCL1
Accession No.	Swiss-Prot#Q07820NCBI Gene ID4170
Cell Localization	Helical
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=3ug/Test ICC=1:50-200 IF=1:50-200

Background

Involved in the regulation of apoptosis versus cell survival, and in the maintenance of viability but not of proliferation. Mediates its effects by interactions with a number of other regulators of apoptosis. Isoform 1 inhibits apoptosis. Isoform 2 promotes apoptosis.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.