

## ZBED2 Antibody FITC Conjugated

Catalog No: #C02668F

Package Size: #C02668F 100ul

## Description

|                       |   |
|-----------------------|---|
| Product Name          | ZBED2 Antibody FITC Conjugated  |
| Host Species          | Rabbit  |
| Clonality             | Polyclonal  |
| Isotype               | IgG   |
| Purification          | Purified by Protein A.  |
| Applications          | ICC IF  |
| Species Reactivity    | Hu  |
| Immunogen Description | KLH conjugated synthetic peptide derived from human ZBED2   |
| Conjugates            | FITC  |
| Target Name           | ZBED2   |
| Other Names           | MGC10796; ZBED2; ZBED2_HUMAN; Zinc finger BED domain containing protein 2; Zinc finger BED domain-containing protein 2; Zinc finger, BED domain containing 2. |
| Concentration         | 1mg ml  |
| Formulation           | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.  |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.   |

## Application Details

ICC=1:50-200 IF=1:50-200

## Background

ZBED2

Note: This product is for in vitro research use only and is not intended for use in humans or animals.