

# RASD1 Antibody HRP Conjugated

Catalog No: #C02898H

Package Size: #C02898H 100ul

## Description

|                       |                                                                                                                  |
|-----------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name          | RASD1 Antibody HRP Conjugated                                                                                    |
| Host Species          | Rabbit                                                                                                           |
| Clonality             | Polyclonal                                                                                                       |
| Isotype               | IgG                                                                                                              |
| Purification          | Purified by Protein A.                                                                                           |
| Applications          | WB IHC-P IHC-F ICC                                                                                               |
| Species Reactivity    | Hu Ms Rt                                                                                                         |
| Immunogen Description | KLH conjugated synthetic peptide aa 30-60 281 derived from human RASD1                                           |
| Conjugates            | HRP                                                                                                              |
| Target Name           | RASD1                                                                                                            |
| Other Names           | AGS1; DEXRAS1; MGC:26290; Dexamethasone-induced Ras-related protein 1; Activator of G-protein signaling 1; RASD1 |
| Accession No.         | Swiss-Prot#Q9Y272NCBI Gene ID51655                                                                               |
| Concentration         | 1mg ml                                                                                                           |
| Formulation           | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                                                 |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.                                  |

## Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200 ICC=1:50-200

## Background

Small GTPase. Negatively regulates the transcription regulation activity of the APBB1 FE65-APP complex via its interaction with APBB1 FE65 (By similarity).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.