

Gastrokine 1 Antibody FITC Conjugated

Catalog No: #C02920F

Package Size: #C02920F 100ul

Description

Product Name	Gastrokine 1 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 130-160 199 derived from human Gastrokine 1
Conjugates	FITC
Target Name	Gastrokine 1
Other Names	FOV; CA11; AMP18; BRICD1; foveolin; Gastrokine-1; 18 kDa antrum mucosa protein; AMP-18; Protein CA11; GKN1; UNQ489 PRO1005
Accession No.	Swiss-Prot#Q9NS71NCBI Gene ID56287
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

Has mitogenic activity and may be involved in maintaining the integrity of the gastric mucosal epithelium.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.