

Uteroglobin Antibody FITC Conjugated

Catalog No: #C02960F

Package Size: #C02960F 100ul

Description

Product Name	Uteroglobin Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-70 96 derived from mouse Uteroglobin
Conjugates	FITC
Target Name	Uteroglobin
Other Names	UG; UGB; Utg; CC10; CC16; CCSP; PCB-BP; Uteroglobin; Clara cell 17 kDa protein; Clara cell phospholipid-binding protein; CCPBP; Clara cells 10 kDa secretory protein; PCB-binding protein; Secretoglobin family 1A member 1; Scgb1a1
Accession No.	Swiss-Prot#Q06318NCBI Gene ID22287
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Binds phosphatidylcholine, phosphatidylinositol, polychlorinated biphenyls (PCB) and weakly progesterone, potent inhibitor of phospholipase A2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.