

CXCL10 Antibody FITC Conjugated

Catalog No: #C03000F

Package Size: #C03000F 100ul

Description

Product Name	CXCL10 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 30-80 98 derived from mouse CXCL10
Conjugates	FITC
Target Name	CXCL10
Other Names	C7; IP1; CRG-2; INP1; IP-1; Ifi1; mob-1; Scyb1; gIP-1; C-X-C motif chemokine 1; 1 kDa interferon gamma-induced protein; Gamma-IP1; Interferon-gamma induced protein CRG-2; Small-inducible cytokine B1; Cxcl1; Crg2
Accession No.	Swiss-Prot#P17515NCBI Gene ID15945
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

In addition to its role as a proinflammatory cytokine, may participate in T-cell effector function and perhaps T-cell development.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.