

KLF17 Antibody FITC Conjugated

Catalog No: #C03450F

Package Size: #C03450F 100ul

Description

Product Name	KLF17 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human ZNF393
Conjugates	FITC
Target Name	KLF17
Other Names	FLJ40160; Krueppel like factor 17; Novel zinc finger protein; RP4 675G8.1; RP4675G8.1; Zfp393; Zinc finger protein 393; ZNF393; KLF17_HUMAN; Krueppel-like factor 17.
Accession No.	NCBI Gene ID128209
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

ICC=1:50-200 IF=1:50-200

Background

KLF17 binds G C-rich sites via its zinc fingers and activates transcription from CACCC-box elements. It may be a germ cell-specific transcription factor that plays important roles in spermatid differentiation and oocyte development.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.