

KCNV2 Antibody FITC Conjugated

Catalog No: #C03629F

Package Size: #C03629F 100ul

Description

Product Name	KCNV2 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 177-210 545 derived from human KCNV2
Conjugates	FITC
Target Name	KCNV2
Other Names	Kv8.2; RCD3B; KV11.1; Potassium voltage-gated channel subfamily V member 2; Voltage-gated potassium channel subunit Kv8.2; KCNV2
Accession No.	Swiss-Prot#Q8TDN2NCBI Gene ID169522
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

Potassium channel subunit. Modulates channel activity by shifting the threshold and the half-maximal activation to more negative values.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.