

ox-LDL Antibody FITC Conjugated

Catalog No: #C03636F

Package Size: #C03636F 100ul

Description

Product Name	ox-LDL Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu
Immunogen Description	Full length plasma protein (Human)
Conjugates	FITC
Target Name	ox-LDL
Other Names	ox-LDL: Cu2SO4 oxidized low density lipoprotein; oxidized low density lipoprotein; Low density lipoprotein; MDA oxidized LDL; MDA oxidized low density lipoprotein.
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Low-density lipoprotein (LDL) is the carrier protein for cholesterol in the blood. LDL binds to its receptor on the capillary walls and thereby mediates the uptake and clearance of cholesterol from the circulation. In atherosclerotic lesions oxidatively modified LDL is found and oxidized LDL is specifically recognized and ingested by macrophages via scavenger receptor A and CD36. Oxidized LDL may be a marker of atherosclerosis but the precise changes in oxidized LDL are not well described. Low-density lipoprotein oxidised with Cu2SO4.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.