

IL-2RB Antibody FITC Conjugated

Catalog No: #C03829F

Package Size: #C03829F 100ul

Description

Product Name	IL-2RB Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 501-551 551 derived from human IL2 Receptor beta
Conjugates	FITC
Target Name	IL2 Receptor beta
Other Names	CD122; IL15RB; P70-75; Interleukin-2 receptor subunit beta; IL-2 receptor subunit beta; IL-2R subunit beta; IL-2RB; High affinity IL-2 receptor subunit beta; p75; IL2RB
Accession No.	Swiss-Prot#P14784NCBI Gene ID3560
Cell Localization	Cytoplasm
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 IF=1:50-200

Background

Receptor for interleukin-2. This beta subunit is involved in receptor mediated endocytosis and transduces the mitogenic signals of IL2.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.