

TLR4 Antibody FITC Conjugated

Catalog No: #C03921F

Package Size: #C03921F 100ul

Description

Product Name	TLR4 Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 155-185 839 derived from human TLR4
Conjugates	FITC
Target Name	TLR4
Other Names	LAC; LPH; LPH1; Lactase-phlorizin hydrolase; Lactase-glycosylceramidase; LCT
Accession No.	Swiss-Prot#P09848NCBI Gene ID3938
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=2ug/Test

Background

LPH splits lactose in the small intestine.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.