

GST tag Antibody HRP Conjugated

Catalog No: #C03959H

Package Size: #C03959H 100ul

Description

Product Name	GST tag Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB
Species Reactivity	Ms GST fusion protein
Crossing Reactivity	GST tag fusion protein
Immunogen Description	Recombinant full length GST tag protein
Conjugates	HRP
Target Name	GST tag
Other Names	glutathione-S-transferase; Glutathione S Transferase; Glutathione S transferase Mu 1; GST 26; GST class mu 1; GST1; GSTM1; GSTM1 1; GSTM1a 1a; GSTM1b 1b; GTH4; GTM1; H B; HB subunit 4; MGC26563; MU 1; MU; SJ26 antigen; SjGST.
Accession No.	NCBI Gene ID2944
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:10000-20000

Background

GST is a 26kDa protein encoded by the parasitic helminth *Schistosoma japonicum* and widely used in the pGEX family of GST plasmid expression vectors as a fusion protein with foreign proteins.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.