

Angiotensin II type 1A receptor Antibody FITC Conjugated

Catalog No: #C03965F

Package Size: #C03965F 100ul

Description

Product Name	Angiotensin II type 1A receptor Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Angiotensin II type 1A receptor : 101-200/359
Conjugates	FITC
Target Name	Angiotensin II type 1A receptor
Other Names	AT1; AG2S; AT1B; AT1R; AT1AR; AT1BR; AT2R1; HAT1R; AGTR1B; Type-1 angiotensin II receptor; Angiotensin II type-1 receptor; AGTR1; AGTR1A; AT2R1B
Accession No.	Swiss-Prot#P30556NCBI Gene ID185
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF 1:100-1:500;

Background

Receptor for angiotensin II. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.