

Angiotensin II type 1A receptor Antibody PE Conjugated

Catalog No: #C03965P

Package Size: #C03965P 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Angiotensin II type 1A receptor Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide derived from human Angiotensin II type 1A receptor: 101-200/359
Conjugates	PE
Target Name	Angiotensin II type 1A receptor
Other Names	AT1; AG2S; AT1B; AT1R; AT1AR; AT1BR; AT2R1; HAT1R; AGTR1B; Type-1 angiotensin II receptor; Angiotensin II type-1 receptor; AGTR1; AGTR1A; AT2R1B
Accession No.	Swiss-Prot#P30556NCBI Gene ID185
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF 1:100-500

Background

Receptor for angiotensin II. Mediates its action by association with G proteins that activate a phosphatidylinositol-calcium second messenger system.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.