

LDL Antibody FITC Conjugated

Catalog No: #C03996F

Package Size: #C03996F 100ul

Description

Product Name	LDL Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms
Immunogen Description	Full length plasma protein (Human)
Conjugates	FITC
Target Name	LDL
Other Names	AD2; LPG; APO-E; LDLQC5; Apolipoprotein E; APOE
Accession No.	Swiss-Prot#P02649NCBI Gene ID348
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Mediates the binding, internalization, and catabolism of lipoprotein particles. It can serve as a ligand for the LDL (apo B E) receptor and for the specific apo-E receptor (chylomicron remnant) of hepatic tissues.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.