

CD162 Antibody HRP Conjugated

Catalog No: #C04029H

Package Size: #C04029H 100ul

Description

Product Name	CD162 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB
Species Reactivity	Hu
Immunogen Description	KLH conjugated synthetic peptide aa 260-300 412 derived from human PSGL-1 CD162
Conjugates	HRP
Target Name	CD162
Other Names	CLA; CD162; PSGL1; PSGL-1; P-selectin glycoprotein ligand 1; Selectin P ligand; SELPLG
Accession No.	Swiss-Prot#Q14242NCBI Gene ID6404
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000

Background

A SLe(x)-type proteoglycan, which through high affinity, calcium-dependent interactions with E-, P- and L-selectins, mediates rapid rolling of leukocytes over vascular surfaces during the initial steps in inflammation. Critical for the initial leukocyte capture.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.