

EPO Antibody HRP Conjugated

Catalog No: #C04046H

Package Size: #C04046H 100ul

Description

Product Name	EPO Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WB IHC-P IHC-F
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 150-193 193 derived from human Erythropoietin
Conjugates	HRP
Target Name	EPO
Other Names	EP; MVCD2; Erythropoietin; EPO
Accession No.	Swiss-Prot#P01588NCBI Gene ID2056
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000 IHC-P=1:50-200 IHC-F=1:50-200

Background

Erythropoietin is the principal hormone involved in the regulation of erythrocyte differentiation and the maintenance of a physiological level of circulating erythrocyte mass.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.