

LAMP2 Antibody PE Conjugated

Catalog No: #C04065P

Package Size: #C04065P 100ul

Description

Product Name	LAMP2 Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	FCM
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 305-355 410 derived from human LAMP2
Conjugates	PE
Target Name	LAMP2
Other Names	LAMPB; CD107b; LAMP-2; LGP110; Lysosome-associated membrane glycoprotein 2; Lysosome-associated membrane protein 2; CD107 antigen-like family member B; LAMP2
Accession No.	Swiss-Prot#P13473NCBI Gene ID3920
Cell Localization	Luminal
Concentration	1mg ml
Formulation	10mM Tris Buffered Saline containing 1% BSA, 50% glycerol and 0.09% sodium azide.
Storage	Store at 4C for 12 months.

Application Details

FCM: 1:20-100

Background

Implicated in tumor cell metastasis. May function in protection of the lysosomal membrane from autodigestion, maintenance of the acidic environment of the lysosome, adhesion when expressed on the cell surface (plasma membrane), and inter- and intracellular signal transduction. Protects cells from the toxic effects of methylating mutagens.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.