

Endothelial lipase Antibody FITC Conjugated

Catalog No: #C04075F

Package Size: #C04075F 100ul

Description

Product Name	Endothelial lipase Antibody FITC Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 310-360 500 derived from human Endothelial lipase
Conjugates	FITC
Target Name	Endothelial lipase
Other Names	EL; EDL; PRO719; Endothelial lipase; Endothelial cell-derived lipase; LIPG; UNQ387 PRO719
Accession No.	Swiss-Prot#Q9Y5X9NCBI Gene ID9388
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

IF=1:50-200

Background

Has phospholipase and triglyceride lipase activities. Hydrolyzes high density lipoproteins (HDL) more efficiently than other lipoproteins. Binds heparin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.