

CD23 Fc EpsilonR II Antibody PE Conjugated

Catalog No: #C04139P

Package Size: #C04139P 100ul

Description

Product Name	CD23 Fc EpsilonR II Antibody PE Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	Flow-Cyt ICC IF
Species Reactivity	Hu Ms Rt
Immunogen Description	KLH conjugated synthetic peptide aa 240-290 321 derived from human CD23
Conjugates	PE
Target Name	CD23 Fc EpsilonR II
Other Names	CD23; FCE2; CD23A; IGEBF; CLEC4J; BLAST-2; Low affinity immunoglobulin epsilon Fc receptor; C-type lectin domain family 4 member J; Fc-epsilon-RII; Immunoglobulin E-binding factor; Lymphocyte IgE receptor; FCER2
Accession No.	Swiss-Prot#P06734NCBI Gene ID2208
Cell Localization	Extracellular, Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

Flow-Cyt=1:50-200 ICC=1:50-200 IF=1:50-200

Background

Low-affinity receptor for immunoglobulin E (IgE) and CR2 CD21. Has essential roles in the regulation of IgE production and in the differentiation of B-cells (it is a B-cell-specific antigen).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.