

Cytochrome P450 1A2 Antibody HRP Conjugated

Catalog No: #C04203H

Package Size: #C04203H 100ul

Description

Product Name	Cytochrome P450 1A2 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WBB B IHC-PB IHC-F
Species Reactivity	HuB MsB RtB
Immunogen Description	KLH conjugated synthetic peptide aa 245-285 513 derived from human Cytochrome P450 1A2
Conjugates	HRP
Target Name	Cytochrome P450 1A2
Other Names	CP12; P3-450; P450(PA); Cytochrome P450 1A2; CYP1A2; Cytochrome P(3)450; Cytochrome P450 4; Cytochrome P450-P3; CYP1A2
Accession No.	Swiss-Prot#P05177NCBI Gene ID1544
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000B B IHC-P=1:50-200B IHC-F=1:50-200B

Background

Cytochromes P450 are a group of heme-thiolate monooxygenases. In liver microsomes, this enzyme is involved in an NADPH-dependent electron transport pathway. It oxidizes a variety of structurally unrelated compounds, including steroids, fatty acids, and xenobiotics. Most active in catalyzing 2-hydroxylation. Caffeine is metabolized primarily by cytochrome CYP1A2 in the liver through an initial N3-demethylation. Also acts in the metabolism of aflatoxin B1 and acetaminophen. Participates in the bioactivation of carcinogenic aromatic and heterocyclic amines. Catalyzes the N-hydroxylation of heterocyclic amines and the O-deethylation of phenacetin.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.