

IL-17C Antibody HRP Conjugated

Catalog No: #C04220H

Package Size: #C04220H 100ul

Description

Product Name	IL-17C Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	WBB B IHC-PB IHC-F
Species Reactivity	HuB MsB RtB B
Immunogen Description	KLH conjugated synthetic peptide aa 90-130 197 derived from human IL-17C
Conjugates	HRP
Target Name	IL-17C
Other Names	CX2; IL-17C; Interleukin-17C; Cytokine CX2; IL17C; UNQ561 PRO1122
Accession No.	Swiss-Prot#Q9P0M4NCBI Gene ID27189
Cell Localization	Secreted
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB=1:500-2000B B IHC-P=1:50-200B IHC-F=1:50-200B

Background

Cytokine that plays a crucial role in innate immunity of the epithelium, including to intestinal bacterial pathogens, in an autocrine manner. Stimulates the production of antibacterial peptides and proinflammatory molecules for host defense by signaling through the NF-kappa-B and MAPK pathways. Acts synergically with IL22 in inducing the expression of antibacterial peptides, including S100A8, S100A9, REG3A and REG3G. Synergy is also observed with TNF and IL1B in inducing DEFB2 from keratinocytes. Depending on the type of insult, may have both protective and pathogenic properties, either by maintaining epithelial homeostasis after an inflammatory challenge or by promoting inflammatory phenotype. Enhanced IL17C IL17RE signaling may also lead to greater susceptibility to autoimmune diseases.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.