

# C5b-9 Antibody HRP Conjugated

Catalog No: #C04256H

Package Size: #C04256H 100ul

## Description

|                       |                                                                                 |
|-----------------------|---------------------------------------------------------------------------------|
| Product Name          | C5b-9 Antibody HRP Conjugated                                                   |
| Host Species          | Rabbit                                                                          |
| Clonality             | Polyclonal                                                                      |
| Isotype               | IgG                                                                             |
| Purification          | Purified by Protein A.                                                          |
| Applications          | IHC-PB IHC-F                                                                    |
| Species Reactivity    | HuB MsB RtB                                                                     |
| Immunogen Description | Purified Human SC5b-9 Complex                                                   |
| Conjugates            | HRP                                                                             |
| Target Name           | C5b-9                                                                           |
| Other Names           | C9D; ARMD15; Complement component C9; C9                                        |
| Accession No.         | Swiss-Prot#P02748NCBI Gene ID735                                                |
| Concentration         | 1mg ml                                                                          |
| Formulation           | 0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.                |
| Storage               | Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |

## Application Details

IHC-P=1:50-200B IHC-F=1:50-200B

## Background

Constituent of the membrane attack complex (MAC) that plays a key role in the innate and adaptive immune response by forming pores in the plasma membrane of target cells. C9 is the pore-forming subunit of the MAC.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.