

CD134 Antibody HRP Conjugated

Catalog No: #C04264H

Package Size: #C04264H 100ul

Description

Product Name	CD134 Antibody HRP Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	Purified by Protein A.
Applications	B IHC-PB IHC-FB ICC
Species Reactivity	HuB MsB RtB B B B
Immunogen Description	KLH conjugated synthetic peptide aa 110-160 277 derived from human TNFRSF4
Conjugates	HRP
Target Name	CD134
Other Names	OX4; ACT35; CD134; IMD16; TXGP1L; Tumor necrosis factor receptor superfamily member 4; ACT35 antigen; OX4L receptor; TAX transcriptionally-activated glycoprotein 1 receptor; TNFRSF4
Accession No.	Swiss-Prot#P43489NCBI Gene ID7293
Cell Localization	Extracellular
Concentration	1mg ml
Formulation	0.01M TBS(pH7.4) with 1% BSA, 0.03% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

B IHC-P=1:50-200B IHC-F=1:50-200B ICC=1:50-200

Background

Receptor for TNFSF4 OX4L GP34. Is a costimulatory molecule implicated in long-term T-cell immunity.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.