

Conjugated PDHA1 (Phospho-Ser293) Rabbit mAb

Catalog No: #C14261

Package Size: #C14261-AF350 100ul #C14261-AF405 100ul #C14261-AF488 100ul #C14261-AF555 100ul #C14261-AF594 100ul
 #C14261-AF647 100ul #C14261-AF680 100ul #C14261-AF750 100ul #C14261-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Conjugated PDHA1 (Phospho-Ser293) Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse Rat
Specificity	Phospho-PDHA1 (S293) Antibody detects endogenous levels of total Phospho-PDHA1 (S293)
Immunogen Description	A synthesized peptide derived from human Phospho-PDHA1 (S293)
Conjugates	AF Dye
Modification	Phospho
Other Names	PDH; PDHA; PDHA1; PDHCE1A; PHE1A; Pyruvate Dehydrogenase (lipoamide) alpha 1; Pyruvate Dehydrogenase E1 alpha;
Accession No.	Uniprot:P08559
Calculated MW	Predicted band size: 43 kDa
SDS-PAGE MW	Observed band size: 40 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4/ Cin dark for 6 months

Application Details

Suggested Dilution::

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.