

Conjugated FGFR3 (Phospho-Tyr724) Rabbit mAb

Catalog No: #C14262

Package Size: #C14262-AF350 100ul #C14262-AF405 100ul #C14262-AF488 100ul #C14262-AF555 100ul #C14262-AF594 100ul
 #C14262-AF647 100ul #C14262-AF680 100ul #C14262-AF750 100ul #C14262-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Conjugated FGFR3 (Phospho-Tyr724) Rabbit mAb
Host Species	Recombinant Rabbit
Clonality	Monoclonal
Purification	Affinity-chromatography
Species Reactivity	Human Mouse
Specificity	Phospho-FGFR3 (Y724) Antibody detects endogenous levels of total Phospho-FGFR3 (Y724)
Immunogen Description	A synthesized peptide derived from human Phospho-FGFR3 (Y724)
Conjugates	AF Dye
Modification	Phospho
Other Names	ACH; CD333; CEK2; FGFR3; Fibroblast growth factor receptor 3; HSFGR3EX; Hydroxyaryl protein kinase; JTK4; Tyrosine kinase JTK 4;
Accession No.	Uniprot:P22607
Calculated MW	Predicted band size: 88 kDa
SDS-PAGE MW	Observed band size: 88 kDa
Formulation	Liquid in PBS, 0.5% BSA, 40% Glycerol, 0.02% Sodium Azide.
Storage	Store at 4/ Cin dark for 6 months

Application Details

Suggested Dilution::

AF350 conjugated: most applications: 1: 50 - 1: 250;

AF405 conjugated: most applications: 1: 50 - 1: 250;

AF488 conjugated: most applications: 1: 50 - 1: 250;

AF555 conjugated: most applications: 1: 50 - 1: 250;

AF594 conjugated: most applications: 1: 50 - 1: 250;

AF647 conjugated: most applications: 1: 50 - 1: 250;

AF680 conjugated: most applications: 1: 50 - 1: 250;

AF750 conjugated: most applications: 1: 50 - 1: 250;

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000;

Note: This product is for in vitro research use only and is not intended for use in humans or animals.